

## References

1. Bunce M, O'Neill CM, Barnardo MCNM, Krausa P, Browning MJ, Morris PJ, Welsh KI. Phototyping: comprehensive DNA typing for HLA-A, B, C, DRB1, DRB3, DRB4, DRB5 & DQB1 by PCR with 144 primer mixes utilising sequence-specific primers (PCR-SSP). *Tissue Antigens* 1995; 46: 355-367.
2. Bodmer JG, Marsh SG, Albert ED, Bodmer WF, Bontrop RE, Charron, Dupont B, Erlich HA, Fauchet R, Bach B, Mayr WR, Parham P, Sasazuki T, Schreuder GM, Strominger JL, Svejgaard A, Terasaki PI. Nomenclature for factors of the HLA System, 1996. *Hum Immunol* 1997, 53: 98-128.
3. Nomenclature for factors of the HLA System. Compiled by Steven G. E. Marsh for the WHO Nomenclature Committee for Factors of the HLA System. <http://www.anthonynolan.com/HIG/nomenc.html>
4. Schaffer M, Olerup O. HLA-AB typing by polymerase-chain reaction with sequence-specific primers: more accurate, less errors, and increased resolution compared to serological typing. *Tissue Antigens*. 2001; 58: 299-307.

# HLA-B\*57 Single Box 1.0 Typing Kit

*In vitro* diagnostics disposal


## Instructions Manual



DESENVOLVIMENTO E PRODUÇÃO  
DE TESTES DE DIAGNÓSTICO

biocant – centro de inovação em biotecnologia  
núcleo 4, lote 3  
3060-197 cantanhede  
portugal

tel + 351 231 410 946  
fax + 351 231 410 947  
e-mail [info@genebox.com](mailto:info@genebox.com)  
[www.genebox.com](http://www.genebox.com)



DESENVOLVIMENTO E PRODUÇÃO  
DE TESTES DE DIAGNÓSTICO

## Table of contents

Presentation .....	2
Product changes and improvements.....	2
Quality Control.....	3
Cell line validation sheet.....	3
HLA-B*57 Single Box 1.0 Typing Kit components.....	4
PCR amplification protocol.....	5
Reagents .....	5
DNA Extraction .....	5
PCR amplification.....	5
PCR cycling parameters.....	6
Gel Electrophoresis protocol .....	7
Preparing 2% agarose gel .....	8
Gel Electrophoresis .....	8
HLA-B*57 Single Box 1.0 plate.....	9
HLA-B*57 Single Box 1.0 plate identification.....	9
Results interpretation sheet.....	10
Results interpretation table.....	10
Troubleshooting Guide .....	11
Precautions and Warnings .....	12
Technical Guide.....	13
Guarantee.....	14
Warranty.....	15
Declaration of Conformity.....	16
Material safety data sheet (MSDS).....	17
References .....	20

## Material Safety Data Sheet (MSDS) (2/3)

### 12. Ecological information

No data available.

### 13. Waste disposal information

Waste disposal in accordance with all applicable regulations (the disposals should be incinerated).

### 14. Transport information

During transportation the temperature could not exceed 25°C.  
Transportation should not exceed 3 days.

### 15. Other information

The above information is based on our current level of knowledge, but does not purport to be all-inclusive and shall be used only as a guide. *geneBOX - R&D Diagnostic Tests* shall not be held liable for any damage resulting from handling or from contact with the above products.

**If your problems persist, do not hesitate to contact our technical support**

**+351 231 410 946**

## Material Safety Data Sheet (MSDS) (2/3)

### 6. Personal protection.

**Hand protection:** Wear appropriate chemically resistant gloves.

**Eye protection:** Chemical safety goggles are recommended.

**Skin protection:** Wear laboratory coat.

### 7. Handling and storage

**Handling:** Avoid substance contact.

**Storage:** Protect from light. Store at temperature indicated on package.

**Package Damage:** reject damaged components.

### 8. Hazards

**Master Mix Components:** may be harmful by **inhalation, ingestion or skin absorption**. May cause eye and skin irritation. Material is irritating to mucous membranes and upper respiratory tract. **Ingestion** of large amounts can cause stomach pains, vomiting or diarrhoea.

### 9. First aid measures

**In case of eye contact:** Immediately flush eyes with large amounts of water for at least 15 minutes. Call a physician.

**In case of skin contact:** Immediately wash skin with soap and large amounts of water. Wash contaminated clothing before re-use.

**In case of ingestion:** Wash out mouth with water provided person is conscious. Call a physician if needed.

**In case of inhalation:** remove to fresh air, if not breathing give artificial respiration. If breathing difficult, give oxygen. Call a physician.

### 10. Fire fighting measures

**Extinguishing media:** Water, carbon dioxide, dry chemical powder or appropriate foam.

**Extinguishing media NOT to use:** None are known.

**Special exposure hazards:** May emit toxic fumes of carbon dioxide, carbon monoxide, nitrogen, phosphorus, hydrogen chloride, and hydrogen gas under fire conditions.

**Special fire-fighting equipment:** When large amounts of substances are released work only with self-contained breathing apparatus and protective clothing to prevent contact with skin and eyes.

### 11. Accidental release measures

**Personal Precautions:** Avoid substance contact. No further requirements.

**Cleaning Method:** Clean up affected area. No further requirements.

## Presentation

Some HLA-B\*57 subtypes are associated with the rate of illness progressions in HIV-1 infected patients, particularly B\*5701 allele in Caucasians and B\*5703 allele in Africans populations. The presence of these alleles is strongly associated with hypersensitivity related with ABACAVIR administration/treatment.

HLA-B\*57 is also related with HCV infection and *psoriasis* susceptibility.

This kit contains typing plates with dried primer mix and PCR Master Mix for rapid screening of HLA-B\*57 gene.

## Product Changes and Improvements

The HLA-B\*57 Single Box specificity and interpretation tables are constantly updated, to include new HLA class I alleles described. This product can also be improved in order to increase the yield of the specific PCR product.

The primers exchanged, added or modified, compared to the previous lot, are detailed in the table below.

Tube	primers	motivation
N/A		

## Quality Control

The following 3 DNA samples from the 13th International Histocompatibility Workshop SSOP Panel were used to verify the kit primers mixes specificity.

Workshop Name	Designation
IHW 09052	DBB
IHW 09273	LADA
IHW 09398	FH18

No false positive or negative amplifications were obtained.

The negative control tube can detect cross-contamination with PCR products.

## Cell line validation sheet

HLA-B*57 low resolution SSP typing kit					
Cell line		Cell Typing			HLA-B*57 Positive well no.
		HLA-A*	HLA-B*	HLA-Cw*	
9273	LADA	0201:8001	0702:5703	0701:0802	All except negative control
9052	DBB	0201	5701	0602	All except negative control
9398	FH18	7401:3601	5301:5703	0401:0701	All except negative control

## Material Safety Data Sheet (MSDS) (1/3)

### geneBOX - R&D Diagnostic Tests™ PCR-SSP Kits

#### geneBOX™ PCR-SSP typing products

This Material Safety Data Sheet (MSDS) applies to all geneBOX - R&D Diagnostic Tests SSP™ typing kits

#### 1. Chemical products and company identification

Date of Issue:	May 2010
Product group:	geneBOX™ PCR-SSP Typing Products
Manufacturer:	geneBOX - R&D Diagnostic Tests, biocant – centro de inovação em biotecnologia núcleo 4, lote 3 3060-197 cantanhede, portugal
tel/fax:	+351 231 410 946/ +351 231 410 947
e-mail:	info@genebox.com

#### 2. Composition and reagents information

Component	Chemical	Common Name
Plate	Deoxyribonucleic acid Cresol Red	Oligonucleotide
PCR Master Mix	Deoxyribonucleotides NH <sub>4</sub> Buffer Magnesium chloride Cresol Red Glycerol	Nucleotides MgCl <sub>2</sub> Glycerine

#### 3. Physic-chemical properties:

Components	Appearance	Colour	Odour
Plate	dried, in plate wells	Red	none
Master Mix	liquid	Pink/red	none

#### 4. Toxicological information

Chemical	Toxicities
Glycerol	LD50= oral 4090 mg/kg (mouse) LD50= oral 12600 mg/kg (rat) LD50= oral 1480 mg/kg (human)

#### 5. Stability and reactivity

**Conditions to avoid:** Heat and moisture.

**Incompatibilities:** Strong oxidizing agents, strong bases.

## Declaration of conformity

**Product Name:** HLA-B\*57 Single Box


**Product Number:** GB.12.07

**Intended use:** HLA-B\*57 low resolution histocompatibility testing.

**Manufacturer:** geneBOX - R&D Diagnostic Tests,  
Biocant – centro de inovação em biotecnologia  
núcleo 4, lote 3  
3060-197 Cantanhede,  
Portugal

We, geneBOX - R&D Diagnostic Tests, hereby declare that this product, to which this declaration of conformity relates, is in conformity with the following standards and other normative documents ISO 9001:2008 and ISO 13485:2003, following the provisions of the 98/79/EC Directive on *in vitro* diagnostic medical devices as transposed into the national laws of the Member States of the European Union.

The technical file of the product is maintained at geneBOX - R&D Diagnostic Tests, Biocant Park, Parque tecnológico de Cantanhede, 3060-197 Cantanhede, Portugal.



Sandra Balseiro  
Technical Director

HLA-B\*57 Single Box 1.0 Typing Kit Components

- HLA-B\*57 Single Box typing plates\*** (88 typing)  
4 plates (22 samples each) (Store at -30 to -15 °C)
- PCR Master Mix (With Taq DNA Polymerase)**  
4 X 80 µl (Store at -30 to -15°C)
- Plate sealers**  
6 sealer caps
- Instructions Manual**  
1 Instructions manual

\* With dried specific primers pairs (1 specific primer pair; 1 negative control and 1 positive control).

PCR Master Mix Components

- Nucleotides**  
Final concentration of each dNTP: 600 µM
- PCR Buffer**  
Final concentration: 3,3x NH<sub>4</sub> Buffer; 2,0 mM MgCl<sub>2</sub> and 0,4 U/µl Taq DNA polymerase, pH 8.3.
- Glycerol**  
Final concentration: 16,6%
- Cresol Red**  
Final concentration: 300µg/ml

aaaaaaa5aaaaaaaaa20a

PCR amplification protocol

Reagents

- DNA Sample (100-200 ng/µl)
- PCR Master Mix
- ddH<sub>2</sub>O (not supplied)

DNA Extraction

For SSP typing highly pure DNA is needed. We recommend isolation of DNA using any extraction kit with CE marking, which guarantees an OD ratio 260/280 higher than 1.6 and a DNA concentration between 100ng – 200ng/µl. Alternatively, the DNA can be extracted using trimethyammoium-bromide salts (DTAB/CTAB) or by salting out, dissolving it in TE Buffer. The same OD and concentration values should be assured.  
DO NOT USE HEPARINISED BLOOD WITH THIS METHOD

PCR Amplification

- Spin briefly the DNA and Master Mix tubes.
- Add:
  - 80 µl of PCR Master Mix,**
  - 160 µl of ddH<sub>2</sub>O**to a 0,7 ml or 1,5 ml tube.
- Vortex the tube vigorously for 15s.
- Load **9 µl** of the mix into each tube of the strip (**1 primers pair**).
- Load **1 µl of DNA sample (conc. 100-200 ng/µl)** in each well (except to the control wells).
- Close the typing plate with a strip cap (sealer cap) and put it in a 96 well thermo cycler.

PCR Cycling Parameters

	Step
	Temperature
	time
	Cycle
	Denaturation
	96 °C
	1 min
	1

aaaaaaa6aaaaaaaaa20a

Precautions and Warnings

PCR amplification allows the amplification of small quantities of sample DNA in an exponential way. However, this is also true for foreign DNA, which can contaminate our PCR method. Consequently, special laboratory practices are necessary in order to avoid false positive amplifications. Bellow is listed Genebox recommendations to circumvent contaminations:

- Work in separated pre-PCR and post-PCR rooms.
- Laboratory workflow must be unidirectional, from pre-PCR to post-PCR area.
- Specific equipment for each working area must be used (sample preparation, amplification and preamplification).
- All equipment used in post-PCR should not leave this area.
- Use dedicated micropipettes, gloves and lab coats in each area.
- Use non talcum powder gloves (since talcum could inhibit the PCR reaction).
- Use filter tips in order to avoid cross contamination.
- Check regularly micropipettes, in order to ensure that they are accurate within 5 % of fixed volume.
- Use different micropipettes depending on the volume we wish to load.
- Check regularly thermocyclers, in order to ensure that they are accurate within 1% of fixed temperature.
- Open and close reagent vials carefully. After use, close vials and store at indicated temperatures.
- Do not use a kit after its expiration date.
- packaging material included within the kit is resistant to the indicated storage conditions. Storage at different conditions can cause breakage of the material, and possible contamination of the kit reagents.
- plastic material included within the kit is resistant under normal conditions of use. Use of plastic material in extreme conditions may cause its breakage, and therefore, the impossibility to use the kit.
- check suitability of DNA quantity and quality before use the kit.

General instructions for laboratory safety:

- do not eat, drink or smoke in laboratory work areas
- wear disposable gloves
- wear clean lab coats and eye protection
- wash hands thoroughly after handling specimens and test reagents
- clean the working area before and after kit handling.
- do not pipette by mouth.

Denaturation  
Annealing  
Extension

96 °C  
70 °C  
72 °C

25 sec  
45 sec  
30 sec

5

Denaturation  
Annealing  
Extension

96 °C  
65 °C  
72 °C

25 sec  
45 sec  
30 sec

21

Denaturation  
Annealing  
Extension

96 °C  
55 °C  
72 °C

25 sec  
1 min  
2 min

4

### Extension

72 °C

aaaaaaa7aaaaaaaaa20a

10 min

1

Keep  
(optional)

4 °C

Infinite

1

7. Keep the plate at 2-8 °C after the PCR have finished.
8. Detect the PCR products with 2% agarose gel electrophoresis.

aaaaaaaa8aaaaaaaaaaaaa20a



Keep different lab coats in pre-PCR and post-PCR rooms
Change protective gloves frequently
Repeat typing with a good quality DNA sample
Blurred bands Degradation of DNA sample Re-extract the DNA sample with fresh material
Repeat typing with a good quality DNA sample
Excess of template DNA Check DNA quality and concentration
Dissolve the DNA sample in ddH2O in order to have the proper concentration
Repeat typing with a good quality DNA sample
Electrophoresis Buffer Problems: wrong buffer or older buffer Use a fresh recommended buffer



1320

Troubleshooting Guide

The control and specific bands are weak. Concentration of DNA sample is too low. Check DNA quality and concentration
Re-extract the DNA sample or try not add water into the PCR Mix
Repeat typing with a good quality DNA sample
DNA polymerase inhibitors in the DNA sample Re-purify the sample DNA
Repeat typing with a good quality DNA sample
Missing internal control bands in one or several lanes. DNA polymerase inhibitors in the DNA sample. Re-purify the sample DNA
Repeat typing with a good quality DNA sample
Dried PCR amplification products Check the plate sealing
Repeat the typing using a PCR MicroMat and/or overlay the PCR reaction mix with mineral oil
False negative of a specific band while the internal control appears normal Degradation of DNA sample Re-extract the DNA sample with fresh material
Repeat typing with a good quality DNA sample
More than two specific alleles are detected/ Ambiguous results Excess of template DNA Check DNA quality and concentration
Dissolve the DNA sample in ddH2O in order to have the proper concentration
Repeat typing with a good quality DNA sample
Contamination with previously amplified PCR products or with other DNA samples during the DNA extraction or PCR preparation steps Clean the working area
Work in separated pre-PCR and post-PCR rooms

1220

Gel Electrophoresis protocol

PREPARING 2% AGAROSE GEL

- 1. Dissolve **4 grams** of electrophoresis **grade agarose** powder in **200 ml** of **1X TAE buffer**.
- 2. Melt the agarose powder completely in a microwave oven.
- 3. Cool the heated agarose gel to ~ **50°C**.
- 4. Add at least 10 µl of **ethidium bromide**<sup>++</sup> (10 mg/ml) or **Sybr Safe**<sup>TM</sup> (100000 x concentrate) to the heated agarose. Stir until it is thoroughly incorporated.
- 5. On a balanced surface, set up a gel plate with **96 wells**.
- 6. Cast a **5mm** thick gel on the plate.
- 7. Allow the gel to settle.


<sup>++</sup>**Caution**, this reagent is a strong mutagenic agent (read carefully its MSDS before using it).

GEL ELECTROPHORESIS

- 1. Submerge the gel in 1X TAE buffer in a gel box.
- 2. Gently remove the caps to avoid splashing of PCR products.
- 3. Load 10 µl into each well on the gel.
- 4. Connect the electric leads and turn on the power supply (115V). Electrophoresis for ~ 20 minutes, or until 2/3 of the lane.
- 5. Transfer the gel onto a UV transilluminator, document the result by photography.
- 6. Use the ***result interpretation sheet*** to interpret results.

aaaaaaaa9aaaaaaaaaaaaa20a

HLA-B\*57 Single Box 1.0 plate



\* Numbers may differ from plate to plate:  
The no 1 could be 4, 7, 10  
The no 2 could be 5, 8, 11  
The no 3 could be 6, 9, 12

HLA-B\*57 Single Box 1.0 Plate Identification

Position
Product
Lines 1, 2 e 3 (except controls)
B*5701-5714 except 5705
CP
Positive control
CN
Negative control

aaaaaaaa10aaaaaaaaaaaaa20a

Results Interpretation sheet

Position HLA-B 57 Specific band Control band**	
Lines 1,2 and 3 (except controls) B*57 positive 380 790 pb	
B*57 negative ---- 790 pb	
PC Positive control 380 ----	
NC Negative control ---- 790 pb	

\*\*Control primer pares match with non-allelic sequences. The internal positive control primer pairs amplify segments of the HLA-DRB1 gene, giving rise to 790 base pair fragment.  
In the presence of the specific band amplification the control band intensity often decreases.  
The PCR reaction is only valid in the presence of control band or, in some cases, in the presence of the specific band.  
In the absence of the control band, please repeat the typing.

1120